

FMNL1 Antibody, FITC conjugated



PACO26799

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in the control of cell motility and survival of macrophages. Plays a role in the regulation of cell morphology and cytoskeletal organization. Required in the cortical actin filament dynamics and cell shape.

Gene ID:

FMNL1

Uniprot

O95466

Synonyms:

Formin-like protein 1 (CLL-associated antigen KW-13) (Leukocyte formin), FMNL1, C17orf1 C17orf1B FMNL FRL1

Immunogen:

Recombinant Human Formin-like protein 1 protein (207-457AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A