

FOSL2 Antibody, HRP conjugated



PACO26814

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Controls osteoclast survival and size. As a dimer with JUN, activates LIF transcription. ctivates CEBPB transcription in PGE2-activated osteoblasts.

Gene ID:

FOSL2

Uniprot

P15408

Synonyms:

Fos-related antigen 2 (FRA-2), FOSL2, FRA2

Immunogen:

Recombinant Human Fos-related antigen 2 protein (1-326AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A