

FOXP1 Antibody, Biotin conjugated



PACO26824

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcriptional repressor. It plays an important role in the specification and differentiation of lung epithelium. Can act with CTBP1 to synergistically repress transcription but CTPBP1 is not essential. Essential transcriptional regulator of B-cell development.

Gene ID:

FOXP1

Uniprot

Q9H334

Synonyms:

Forkhead box protein P1 (Mac-1-regulated forkhead) (MFH), FOXP1

Immunogen:

Recombinant Human Forkhead box protein P1 protein (1-114AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A