

FRS3 Antibody, Biotin conjugated



PACO26836

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Adapter protein that links FGF and NGF receptors to downstream signaling pathways. Involved in the activation of MAP kinases. Down-regulates ERK2 signaling by interfering with the phosphorylation and nuclear translocation of ERK2.

Gene ID:

FRS3

Uniprot

O43559

Synonyms:

Fibroblast growth factor receptor substrate 3 (FGFR substrate 3) (FGFR-signaling adaptor SNT2) (Suc1-associated neurotrophic factor target 2) (SNT-2), FRS3

Immunogen:

Recombinant Human Fibroblast growth factor receptor substrate 3 protein (101-400AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A