

FSCN2 Antibody



PACO26837

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Acts as an actin bundling protein. May play a pivotal role in photoreceptor cell-specific events, such as disk morphogenesis.

Gene ID:

FSCN2

Uniprot

O14926

Synonyms:

Fascin-2 (Retinal fascin), FSCN2

Immunogen:

Recombinant Human Fascin-2 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A