

# FTSJ1 Antibody, FITC conjugated



PACO26843

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FTSJ1

**Uniprot**

Q9UET6

**Synonyms:**

Putative tRNA (cytidine(32)/guanosine(34)-2'-O)-methyltransferase (EC 2.1.1.205) (2'-O-ribose RNA methyltransferase TRM7 homolog) (Protein ftsJ homolog 1), FTSJ1

**Immunogen:**

Recombinant Human Putative tRNA (cytidine(32)/guanosine(34)-2'-O)-methyltransferase protein (1-329AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A