

FXYP7 Antibody



PAC026857

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

integral component of plasma membrane, plasma membrane, sodium channel regulator activity, ion transmembrane transport, regulation of sodium ion transmembrane transporter activity, transmembrane transport.

Gene ID:

FXYP7

Uniprot

P58549

Synonyms:

FXYP domain-containing ion transport regulator 7, FXYP7


Immunogen:

Recombinant Human FXYP domain-containing ion transport regulator 7 protein (47-80AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human heart tissue using PACO26857 at dilution of 1:100.



Immunofluorescent analysis of U251 cells using PACO26857 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).