
Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the production of GABA.

Gene ID:

Gad2

Uniprot

P48320

Synonyms:

Glutamate decarboxylase 2 (EC 4.1.1.15) (65 kDa glutamic acid, decarboxylase) (GAD-65) (Glutamate decarboxylase 65 kDa isoform), Gad2, Gad65

Immunogen:

Recombinant Mouse Glutamate decarboxylase 2 protein (1-585AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A