

PACO26889

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

In T-cells, functions as a regulator of p38 MAPKs by inhibiting p38 phosphorylation and activity. Might affect PCNA interaction with some CDK (cell division protein kinase) complexes; stimulates DNA excision repair in vitro and inhibits entry of cells into S phase.

Gene ID:

GADD45A

Uniprot

P24522

Synonyms:

Growth arrest and DNA damage-inducible protein GADD45 alpha (DNA damage-inducible transcript 1 protein) (DDIT-1), GADD45A, DDIT1 GADD45

Immunogen:

Recombinant Human Growth arrest and DNA damage-inducible protein GADD45 α protein (1-165AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A