

# GADD45A Antibody, Biotin conjugated



PACO26892

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

In T-cells, functions as a regulator of p38 MAPKs by inhibiting p38 phosphorylation and activity. Might affect PCNA interaction with some CDK (cell division protein kinase) complexes; stimulates DNA excision repair in vitro and inhibits entry of cells into S phase.

**Gene ID:**

GADD45A

**Uniprot**

P24522

**Synonyms:**

Growth arrest and DNA damage-inducible protein GADD45 alpha (DNA damage-inducible transcript 1 protein) (DDIT-1), GADD45A, DDIT1 GADD45

**Immunogen:**

Recombinant Human Growth arrest and DNA damage-inducible protein GADD45 &alpha; protein (1-165AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A