

GADD45G Antibody, HRP conjugated



PACO26898

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the regulation of growth and apoptosis. Mediates activation of stress-responsive MTK1/MEKK4 MAPKKK.

Gene ID:

GADD45G

Uniprot

O95257

Synonyms:

Growth arrest and DNA damage-inducible protein GADD45 gamma (Cytokine-responsive protein CR6) (DNA damage-inducible transcript 2 protein) (DDIT-2), GADD45G, CR6 DDIT2

Immunogen:

Recombinant Human Growth arrest and DNA damage-inducible protein GADD45 γ protein (1-159AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A