

GAS2 Antibody, Biotin conjugated



PACO26928

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play a role in apoptosis by acting as a cell death substrate for caspases. Is cleaved during apoptosis and the cleaved form induces dramatic rearrangements of the actin cytoskeleton and potent changes in the shape of the affected cells. May be involved in the membrane ruffling process.

Gene ID:

GAS2

Uniprot

O43903

Synonyms:

Growth arrest-specific protein 2 (GAS-2), GAS2

Immunogen:

Recombinant Human Growth arrest-specific protein 2 protein (1-141AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A