

GCNT2 Antibody, Biotin conjugated



PACO26932

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Branching enzyme that converts linear into branched poly-N-acetyllactosaminoglycans. Introduces the blood group I antigen during embryonic development. It is closely associated with the development and maturation of erythroid cells. Isoform C: Determines the expression of the blood group I antigen in erythrocytes.

Gene ID:

GCNT2

Uniprot

Q8N0V5

Synonyms:

N-acetyllactosaminide beta-1,6-N-acetylglucosaminyl-transferase (N-acetylglucosaminyltransferase) (I-branching enzyme) (IGNT), GCNT2, GCNT5 II NACGT1

Immunogen:

Recombinant Human N-acetyllactosaminide beta-1,6-N-acetylglucosaminyl-transferase protein (201-402AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A