

# GGTLC2 Antibody, FITC conjugated



PACO26967

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

anchored component of external side of plasma membrane, extracellular exosome, gamma-glutamyltransferase activity, glutathione metabolic process, leukotriene biosynthetic process.

**Gene ID:**

GGTLC2

**Uniprot**

Q14390

**Synonyms:**

Glutathione hydrolase light chain 2 (Gamma-glutamyltransferase light chain 2) (Gamma-glutamyltransferase-like protein 4), GGTLC2, GGTL4

**Immunogen:**

Recombinant Human Glutathione hydrolase light chain 2 protein (1-218AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A