

GJB6 Antibody, FITC conjugated



PACO26999

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

One gap junction consists of a cluster of closely packed pairs of transmembrane channels, the connexons, through which materials of low MW diffuse from one cell to a neighboring cell.

Gene ID:

GJB6

Uniprot

O95452

Synonyms:

Gap junction beta-6 protein (Connexin-30) (Cx30), GJB6

Immunogen:

Recombinant Human Gap junction beta-6 protein (216-261AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A