

# GK2 Antibody, HRP conjugated



PACO27006

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Key enzyme in the regulation of glycerol uptake and metabolism.

**Gene ID:**

GK2

**Uniprot**

Q14410

**Synonyms:**

Glycerol kinase 2 (GK 2) (Glycerokinase 2) (EC 2.7.1.30) (ATP: glycerol 3-phosphotransferase 2) (Glycerol kinase, testis specific 2), GK2, GKP2 GKTA

**Immunogen:**

Recombinant Human Glycerol kinase 2 protein (2-553AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A