

PACO27053

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

enzyme activator activity, protein kinase inhibitor activity, negative regulation of cyclin-dependent protein kinase activity, positive regulation of catalytic activity, protein phosphorylation.

**Gene ID:**

GMFG

**Uniprot**

O60234

**Synonyms:**

Glia maturation factor gamma (GMF-gamma), GMFG

**Immunogen:**

Recombinant Human Glia maturation factor &gamma; protein (1-142AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A