

GMFG Antibody, HRP conjugated



PACO27054

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

enzyme activator activity, protein kinase inhibitor activity, negative regulation of cyclin-dependent protein kinase activity, positive regulation of catalytic activity, protein phosphorylation.

Gene ID:

GMFG

Uniprot

O60234

Synonyms:

Glia maturation factor gamma (GMF-gamma), GMFG

Immunogen:

Recombinant Human Glia maturation factor γ protein (1-142AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A