

GNA14 Antibody, FITC conjugated



PACO27063

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling systems.

Gene ID:

GNA14

Uniprot

O95837

Synonyms:

Guanine nucleotide-binding protein subunit alpha-14 (G alpha-14) (G-protein subunit alpha-14), GNA14

Immunogen:

Recombinant Human Guanine nucleotide-binding protein subunit alpha-14 protein (8-303AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A