

GNL1 Antibody, HRP conjugated



PACO27102

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Possible regulatory or functional link with the histocompatibility cluster.

Gene ID:

GNL1

Uniprot

P36915

Synonyms:

Guanine nucleotide-binding protein-like 1 (GTP-binding protein HSR1), GNL1, HSR1

Immunogen:

Recombinant Human Guanine nucleotide-binding protein-like 1 protein (251-500AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A