

Gpc4 Antibody, HRP conjugated



PACO27130

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cell surface proteoglycan that bears heparan sulfate. May be involved in the development of kidney tubules and of the central nervous system.

Gene ID:

Gpc4

Uniprot

P51655

Synonyms:

Glypican-4 (K-glypican) [Cleaved into: Secreted glypican-4], Gpc4

Immunogen:

Recombinant Mouse Glypican-4 protein (530-557AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A