

# GPC4 Antibody, FITC conjugated



PACO27135

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Cell surface proteoglycan that bears heparan sulfate. May be involved in the development of kidney tubules and of the central nervous system.

**Gene ID:**

GPC4

**Uniprot**

O75487

**Synonyms:**

Glypican-4 (K-glypican) [Cleaved into: Secreted glypican-4], GPC4

**Immunogen:**

Recombinant Human Glypican-4 protein (19-210AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A