

GPD2 Antibody, FITC conjugated



PACO27147

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

mitochondrial inner membrane, calcium ion binding, glycerol-3-phosphate dehydrogenase activity, cellular lipid metabolic process, small molecule metabolic process.

Gene ID:

GPD2

Uniprot

P43304

Synonyms:

Glycerol-3-phosphate dehydrogenase, mitochondrial (GPD-M) (GPDH-M) (EC 1.1.5.3) (mtGPD), GPD2

Immunogen:

Recombinant Human Glycerol-3-phosphate dehydrogenase, mitochondrial protein (43-727AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A