

PPP2CB Antibody, HRP conjugated



PACO27174

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

PP2A can modulate the activity of phosphorylase B kinase casein kinase 2, mitogen-stimulated S6 kinase, and MAP-2 kinase.

Gene ID:

PPP2CB

Uniprot

P62714

Synonyms:

Serine/threonine-protein phosphatase 2A catalytic subunit beta isoform (PP2A-beta) (EC 3.1.3.16), PPP2CB

Immunogen:

Recombinant Human Serine/threonine-protein phosphatase 2A catalytic subunit β isoform protein (1-309AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A