

GPX7 Antibody, Biotin conjugated



PACO27184

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

It protects esophageal epithelia from hydrogen peroxide-induced oxidative stress. It suppresses acid, c bile acid, induced reactive oxygen species (ROS) and protects against oxidative DNA damage and double-strand breaks.

Gene ID:

GPX7

Uniprot

Q96SL4

Synonyms:

Glutathione peroxidase 7 (GPx-7) (GSHPx-7) (EC 1.11.1.9) (CL683), GPX7, GPX6

Immunogen:

Recombinant Human Glutathione peroxidase 7 protein (28-187AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A