

GSDMD Antibody, FITC conjugated



PACO27207

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

nucleoplasm.

Gene ID:

GSDMD

Uniprot

P57764

Synonyms:

Gasdermin-D (Gasdermin domain-containing protein 1) [Cleaved into: Gasdermin-D, N-terminal (GSDMD-NT); Gasdermin-D, C-terminal (GSDMD-CT)], GSDMD, DFNA5L, GSDMDC1

Immunogen:

Recombinant Human Gasdermin-D protein (1-256AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A