

GSTA4 Antibody, HRP conjugated



PACO27222

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. This isozyme has a high catalytic efficiency with 4-hydroxyalkenals such as 4-hydroxynonenal (4-HNE).

Gene ID:

GSTA4

Uniprot

O15217

Synonyms:

Glutathione S-transferase A4 (EC 2.5.1.18) (GST class-alpha member 4) (Glutathione S-transferase A4-4), GSTA4

Immunogen:

Recombinant Human Glutathione S-transferase A4 protein (1-222AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A