

GSTK1 Antibody, HRP conjugated



PACO27226

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Significant glutathione conjugating activity is found only with the model substrate, 1-chloro-2,4-dinitrobenzene (CDNB).

Gene ID:

GSTK1

Uniprot

Q9Y2Q3

Synonyms:

Glutathione S-transferase κ 1 (EC 2.5.1.18) (GST 13-13) (GST class-kappa) (GSTK1-1) (hGSTK1) (Glutathione S-transferase subunit 13), GSTK1

Immunogen:

Recombinant Human Glutathione S-transferase κ 1 protein (2-226AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A