

# GTF2F2 Antibody, FITC conjugated



PACO27255

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

TFIIF is a general transcription initiation factor that binds to RNA polymerase II and helps to recruit it to the initiation complex in collaboration with TFIIB. It promotes transcription elongation. This subunit shows ATP-dependent DNA-helicase activity.

**Gene ID:**

GTF2F2

**Uniprot**

P13984

**Synonyms:**

General transcription factor IIF subunit 2 (EC 3.6.4.12) (ATP-dependent helicase GTF2F2) (General transcription factor IIF 30 kDa subunit) (Transcription initiation factor IIF subunit beta) (TFIIF-beta) (Transcription initiation factor RAP30), GTF2F2, RAP30

**Immunogen:**

Recombinant Human General transcription factor IIF subunit 2 protein (2-249AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A