

GUCA1A Antibody, HRP conjugated



PACO27262

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Stimulates guanylyl cyclase 1 (GC1) when free calcium ions concentration is low and inhibits GC1 when free calcium ions concentration is elevated. This Ca(2+)-sensitive regulation of GC is a key event in recovery of the dark state of rod photoreceptors following light exposure.

Gene ID:

GUCA1A

Uniprot

P43080

Synonyms:

Guanylyl cyclase-activating protein 1 (GCAP 1) (Guanylate cyclase activator 1A), GUCA1A, C6orf131 GCAP GCAP1 GUCA1

Immunogen:

Recombinant Human Guanylyl cyclase-activating protein 1 protein (2-201AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A