

PACO27265

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Endogenous activator of intestinal guanylate cyclase. It stimulates this enzyme through the same receptor binding region as the heat-stable enterotoxins. May be a potent physiological regulator of intestinal fluid and electrolyte transport. May be an autocrine/paracrine regulator of intestinal salt and water transport.

**Gene ID:**

Guca2b

**Uniprot**

O09051

**Synonyms:**

Guanylate cyclase activator 2B [Cleaved into: Uroguanylin (UGN)], Guca2b

**Immunogen:**

Recombinant Mouse Guanylate cyclase activator 2B protein (22-106AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A