

Guca2b Antibody, Biotin conjugated



PACO27268

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Endogenous activator of intestinal guanylate cyclase. It stimulates this enzyme through the same receptor binding region as the heat-stable enterotoxins. May be a potent physiological regulator of intestinal fluid and electrolyte transport. May be an autocrine/paracrine regulator of intestinal salt and water transport.

Gene ID:

Guca2b

Uniprot

O09051

Synonyms:

Guanylate cyclase activator 2B [Cleaved into: Uroguanylin (UGN)], Guca2b

Immunogen:

Recombinant Mouse Guanylate cyclase activator 2B protein (22-106AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A