

# Hba Antibody



PACO27313

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse, Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:1000-1:5000

**Protein Background:****Gene ID:**

Hba

**Uniprot**

P01942

**Synonyms:**

Hemoglobin subunit alpha (Alpha-globin) (Hemoglobin alpha chain), Hba, Hba-a1

**Immunogen:**


Recombinant Mouse Hemoglobin subunit &alpha; protein (2-142AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



### Western blot

All lanes: Hba antibody at 2 $\mu$ g/ml

Lane 1: HepG2 whole cell lysate

Lane 2: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 16 kDa

Observed band size: 16 kDa

.