

HBE1 Antibody, HRP conjugated



PACO27326

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

The epsilon chain is a beta-type chain of early mammalian embryonic hemoglobin.

Gene ID:

HBE1

Uniprot

P02100

Synonyms:

Hemoglobin subunit epsilon (Epsilon-globin) (Hemoglobin epsilon chain), HBE1, HBE

Immunogen:

Recombinant Human Hemoglobin subunit epsilon protein (2-147AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A