

HBG1 Antibody, HRP conjugated



PACO27330

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

gamma chains make up the fetal hemoglobin F, in combination with alpha chains.

Gene ID:

HBG1

Uniprot

P69891

Synonyms:

Hemoglobin subunit gamma-1 (Gamma-1-globin) (Hb F Agamma) (Hemoglobin gamma-1 chain) (Hemoglobin gamma-A chain), HBG1

Immunogen:

Recombinant Human Hemoglobin subunit gamma-1 protein (2-147AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A