

# HBG1 Antibody, FITC conjugated



PACO27331

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

gamma chains make up the fetal hemoglobin F, in combination with alpha chains.

**Gene ID:**

HBG1

**Uniprot**

P69891

**Synonyms:**

Hemoglobin subunit gamma-1 (Gamma-1-globin) (Hb F Agamma) (Hemoglobin gamma-1 chain) (Hemoglobin gamma-A chain), HBG1

**Immunogen:**

Recombinant Human Hemoglobin subunit gamma-1 protein (2-147AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A