

DAP Antibody, FITC conjugated



PACO27359

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in mediating interferon-gamma-induced cell death.

Gene ID:

DAP

Uniprot

P51397

Synonyms:

Death-associated protein 1 (DAP-1), DAP, DAP1

Immunogen:

Recombinant Human Death-associated protein 1 protein (2-102AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A