

HMOX1 Antibody, HRP conjugated



PACO27398

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Heme oxygenase cleaves the heme ring at the alpha methene bridge to form biliverdin. Biliverdin is subsequently converted to bilirubin by biliverdin reductase. Under physiological conditions, the activity of heme oxygenase is highest in the spleen, where senescent erythrocytes are sequestered and destroyed.

Gene ID:

HMOX1

Uniprot

P09601

Synonyms:

Heme oxygenase 1 (HO-1) (EC 1.14.14.18), HMOX1, HO HO1

Immunogen:

Recombinant Human Heme oxygenase 1 protein (3-288AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A