

HNRNPA2B1 Antibody, FITC conjugated



PACO27407

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved with pre-mRNA processing. Forms complexes (ribonucleosomes) with at least 20 other different hnRNP and heterogeneous nuclear RNA in the nucleus.

Gene ID:

HNRNPA2B1

Uniprot

P22626

Synonyms:

Heterogeneous nuclear ribonucleoproteins A2/B1 (hnRNP A2/B1), HNRNPA2B1, HNRPA2B1

Immunogen:

Recombinant Human Heterogeneous nuclear ribonucleoproteins A2/B1 protein (1-249AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A