

HOXC11 Antibody, FITC conjugated



PACO27431

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. Binds to a promoter element of the lactase-phlorizin hydrolase gene.

Gene ID:

HOXC11

Uniprot

O43248

Synonyms:

Homeobox protein Hox-C11 (Homeobox protein Hox-3H), HOXC11, HOX3H

Immunogen:

Recombinant Human Homeobox protein Hox-C11 protein (10-304AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A