

HOXC9 Antibody, Biotin conjugated



PACO27436

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Gene ID:

HOXC9

Uniprot

P31274

Synonyms:

Homeobox protein Hox-C9 (Homeobox protein Hox-3B), HOXC9, HOX3B

Immunogen:

Recombinant Human Homeobox protein Hox-C9 protein (18-260AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A