

# IDH3G Antibody, HRP conjugated



PACO27490

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

IDH3G

**Uniprot**

P51553

**Synonyms:**

Isocitrate dehydrogenase [NAD] subunit gamma, mitochondrial (Isocitric dehydrogenase subunit gamma) (NAD(+)-specific ICDH subunit gamma), IDH3G

**Immunogen:**

Recombinant Human Isocitrate dehydrogenase [NAD] subunit gamma, mitochondrial protein (40-393AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A