

IDO1 Antibody, Biotin conjugated



PACO27496

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

IDO1

Uniprot

P14902

Synonyms:

Indoleamine 2,3-dioxygenase 1 (IDO-1) (EC 1.13.11.52) (Indoleamine-pyrrole 2,3-dioxygenase), IDO1, IDO INDO

Immunogen:

Recombinant Human Indoleamine 2,3-dioxygenase 1 protein (1-403AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A