

IDUA Antibody, HRP conjugated



PACO27502

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

IDUA

Uniprot

P35475

Synonyms:

Alpha-L-iduronidase (EC 3.2.1.76), IDUA

Immunogen:

Recombinant Human Alpha-L-iduronidase protein (28-653AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A