

IFI44 Antibody, Biotin conjugated



PACO27508

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This protein aggregates to form microtubular structures.

Gene ID:

IFI44

Uniprot

Q8TCB0

Synonyms:

Interferon-induced protein 44 (p44) (Microtubule-associated protein 44), IFI44, MTAP44

Immunogen:

Recombinant Human Interferon-induced protein 44 protein (101-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A