

IFNA14 Antibody, Biotin conjugated



PACO27532

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.

Gene ID:

IFNA14

Uniprot

P01570

Synonyms:

Interferon alpha-14 (IFN-alpha-14) (Interferon alpha-H) (LeIF H) (Interferon lambda-2-H), IFNA14

Immunogen:

Recombinant Human Interferon alpha-14 protein (24-189AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A