

PGD Antibody, Biotin conjugated



PACO27572

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the oxidative decarboxylation of 6-phosphogluconate to ribulose 5-phosphate and CO₂, with concomitant reduction of NADP to NADPH.

Gene ID:

PGD

Uniprot

P52209

Synonyms:

6-phosphogluconate dehydrogenase, decarboxylating (EC 1.1.1.44), PGD, PGDH

Immunogen:

Recombinant Human 6-phosphogluconate dehydrogenase, decarboxylating protein (4-483AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A