

IL10 Antibody, Biotin conjugated



PACO27584

Product Information

Size:

50ug

Reactivity:

Macaca mulatta

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Inhibits the synthesis of a number of cytokines, including IFN-gamma, IL-2, IL-3, TNF and GM-CSF produced by activated macrophages and by helper T-cells.

Gene ID:

IL10

Uniprot

P51496

Synonyms:

Interleukin-10 (IL-10) (Cytokine synthesis inhibitory factor) (CSIF), IL10

Immunogen:

Recombinant Macaca mulatta Interleukin-10 protein (19-178AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A