

IL1B Antibody




PACO27597

Product Information

Size: 50ug	Protein Background:
Reactivity: Sheep	Gene ID: IL1B
Source: Rabbit	Uniprot P21621
Isotype: IgG	Synonyms: Interleukin-1 beta (IL-1 beta), IL1B
Applications: ELISA, WB	Immunogen: Recombinant Sheep Interleukin-1 β protein (114-266AA).
Recommended dilutions: ELISA:1:2000-1:10000, WB:1:1000-1:5000	Storage: Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: IL1B antibody at 2 μ g/ml

Lane 1: Recombinant Interleukin-1 beta protein 0.05 μ g

Lane 2: Recombinant Interleukin-1 beta protein 0.5 μ g

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 31 kDa

Observed band size: 31 kDa