

# IL1B Antibody



PACO27601

---

## Product Information

**Size:**

50ug

**Reactivity:**

Sheep

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:1000-1:5000

**Protein Background:**

Produced by activated macrophages, IL-1 stimulates thymocyte proliferation by inducing IL-2 release, B-cell maturation and proliferation, and fibroblast growth factor activity. IL-1 proteins are involved in the inflammat.

**Gene ID:**

IL1B

**Uniprot**

P21621

**Synonyms:**

Interleukin-1 beta (IL-1 beta), IL1B

**Immunogen:**


Recombinant Sheep Interleukin-1 &beta; protein (114-266AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



### Western blot

All lanes: IL1B antibody at 2 $\mu$ g/ml + Recombinant Interleukin-1 beta protein 0.01 $\mu$ g

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 31 kDa

Observed band size: 31 kDa